


THE WEANING TABLE* - Plan view

A stable low wooden table with a steel plate underneath to avoid tipping over, when the child pulls up on it.



3/4" Baltic top (.75") and apron with poplar legs.
Legs are 2x2" and 12.75" tall

*This table is first used with the weaning chair and is for the weaning age child who starting to sit at mealtime when new foods are being introduced. It is important to give attention to the presentation of the meal and set the table properly with a table cloth or place mat, bib, napkin, plate, cup, jar, flowers etc. A tablecloth should be embroidered with a space for the plate, glass, fork, and spoon.